

# LabCorp At Home COVID-19 Test Home Collection Kit



**1**  
**Read through the instructions completely before collecting your sample.** To see an instructional video or report a missing item in the kit, visit [www.LabCorp.com/AtHomeKit](http://www.LabCorp.com/AtHomeKit).



**2**  
**Your specimen must be shipped Monday through Friday and on the same day that you collect it.** Do not collect your specimen until you confirm pickup times of the drop box or have confirmed your scheduled pickup.  
Visit <https://bit.ly/fedexdropbox> to view FedEx drop box locations and pickup schedules. For more information, refer to the back of this instruction card.




**3**  
**Wash and dry hands before opening the kit. Open your kit and place all the contents on a clean, dry surface.**




**4**  
**Take the swab out of its package.** Do not touch the tip of the swab with your hands.




**5**  
**Screw off the top of the collection tube.** Hold swab in one hand and collection tube in the other, being careful not to spill the liquid. Do not drink the liquid.



**6**  
**Insert the tip of the swab into one nostril.** The swab does not need to be inserted far—insert just until the tip of the swab is no longer visible. Rotate the swab in a circle around the entire inside edge of your nostril at least 3 times.



**7**  
**Take the swab out of your nostril. Using the same end of the swab, repeat step 6 in your other nostril.**



**8**  
**Remove the swab from your second nostril and place in the collection tube.** The end of the swab that went into your nose should be placed into the tube first so that it sits down in the liquid. Screw the top of the collection tube back on.



**9**  
**Wash and dry hands thoroughly again. Insert collection tube into the biohazard specimen bag.** Seal the biohazard specimen bag. Fold the specimen bag in half.



**10**  
**Fill out the Specimen Confirmation Form.** Provide the date and time of collection, and complete the Personal Information section.



**11**  
**Place the specimen bag and Specimen Confirmation Form into the white shipping box and close the lid.** Place the shipping box into the FedEx return pack. Remove the adhesive cover strip and seal the FedEx return pack closed.



**12**  
**Return the postage-paid, pre-addressed FedEx return pack to LabCorp.** You can deliver it to a FedEx drop box or schedule a pickup. It's important to ship your specimen Monday through Friday and on the same day you collect it.



**13**  
**Results will be available in your LabCorp Patient™ account 1 to 2 days after your specimen arrives at the lab.** If you don't have an account, visit [patient.LabCorp.com](http://patient.LabCorp.com) to create one. Results will also be delivered to the ordering healthcare provider.

- This home collection kit has not been FDA cleared or approved.
- This home collection kit has been authorized by FDA under an EUA. This home collection kit has been authorized only for the home collection and maintenance of nasal swab specimens as an aid in detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens.
- This home collection kit in combination with the authorized test is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostic tests for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner.